



- The AltoSa-SFT Balloon is especially designed for large lumen dilatation and tacking of endoprosthesis.
- The short deflation time makes the AltoSa-SFT safe and easy to use.
- **HD** (PU) Material makes the AltoSa-SFT Tacking balloon highly puncture resistant and reliable for multiple dilatations.
- Increased visibility leading to more accurate positioning.

Technical Information:

Balloon Material:	HD PU Material
Balloon Compliance:	Compliant
Shaft Diameter:	10.5 F
Marker Distance & Balloon Length:	40mm & max 46mm
Recommended max Diameter:	46 mm at 70ml
Rated Burst:	100 ml / 50mm
Shaft Length:	65 cm / 100 cm
Recommended Guide Wire	0.035" (0.89 mm)
Recommended Sheath Introducer:	12F

* With this Balloon AndraTec recommends an Extra stiff guide wire such as the Lokum Lunderquist.

ALTOSA-SFT Tacking Catheter Order Codes

Balloon Length Diameter / CSI / F	40mm 100cm shaft	40mm 65cm shaft
46 mm / 12F	-	AL-46-65



Distributed By:



AndraTec GmbH
 Simmernerstr. 70
 D-56075 Koblenz
 Tel: +49 261 29675 172
 Fax: +49 261 29675 173
 www.andratec.com