


AndraTec *Advanced Interventional Technologies*

The **AndraTec** *Exeter* Basket Range features all the requirements that are needed for foreign body extraction. The broad availability makes it easy to choose the right device for the physicians needs. All **AndraTec** Baskets are made from highly kink-resistant Nitinol material that offer a superior strength. Furthermore, **AndraTec** has optimized the helical / Straight shape enabling to easy manipulate and retrieve foreign bodies.



Exeter Retrieval Basket Laser (co-axial) Order Information				
Basket Shape	French Size (F)	Diameter	Length	Order Code
Straight 4 Struts	3F	12 mm	90 cm	COL-30F-1290
Straight 4 Struts	3F	16 mm	90 cm	COL-30F-1690
Straight 4 Struts	3F	12 mm	120 cm	COL-30F-1220
Straight 4 Struts	3F	16 mm	120 cm	COL-30F-1620

Exeter Retrieval Basket Tipless Straight Order Information				
Basket Shape	French Size (F)	Diameter	Length	Order Code
Straight 4 Struts	1.3F	8 mm	120 cm	BTL-13F-0820
Straight 4 Struts	1.3F	12 mm	120 cm	BTL-13F-1220
Straight 4 Struts	1.9F	8 mm	120 cm	BTL-19F-0820
Straight 4 Struts	1.9F	12 mm	120 cm	BTL-19F-1220
Straight 4 Struts	1.9F	16 mm	120 cm	BTL-19F-1620
Straight 4 Struts	2.2F	12 mm	120 cm	BTL-22F-1220
Straight 4 Struts	2.5F	12 mm	120 cm	BTL-25F-1220
Straight 4 Struts	2.5F	16 mm	120 cm	BTL-25F-1620
Straight 4 Struts	3F	12 mm	90 cm	BTL-30F-1290
Straight 4 Struts	3F	16 mm	90 cm	BTL-30F-1690

Exeter Basket Tipless Straight Flat Wire Order Information				
Basket Shape	French Size (F)	Diameter	Length	Order Code
Straight 4 Struts	1.9F	12 mm	120 cm	BTF-19F-1220

Exeter Retrieval Basket With Tip Straight Order Information				
Basket Shape	French Size (F)	Diameter	Length	Order Code
Helical 4 Struts	2.5F	16 mm	90 cm	BWT-25F-1690
Helical 4 Struts	2.5F	16 mm	120 cm	BWT-25F-1620
Helical 4 Struts	3.0F	16 mm	90 cm	BWT-30F-1690
Helical 4 Struts	3.0F	16 mm	120 cm	BWT-30F-1620
Helical 4 Struts	4.0F	16 mm	90 cm	BWT-40F-1690
Helical 4 Struts	4.0F	16 mm	120 cm	BWT-40F-1620



Distributed By:



AndraTec GmbH
 Simmerstr. 70
 D-56075 Koblenz
 Tel: +49 261 29675 172
 Fax: +49 261 29675 173
 www.andratec.com