



Lokum
Lunderquist Guide Wire

- The LOKUM Lunderquist Guide Wire is especially designed for extra strong support and high pushability in difficult procedures.
- The high performance PTFE coating combined with extra stiff stainless steel LOKUM Lunderquist is an excellent choice to treat difficult vascular and non-vascular lesions.
- The Lokum Lunderquist Type guide wire is intended to use for complex diagnostic and interventional procedures where increased stability, and low surface friction of the guidewire is needed.

Technical Information:

| | |
|----------------------|----------------------------------|
| Guide Wire Material: | Extra Tempered Stainless Steel |
| Coating: | PTFE |
| Diameter: | 0.035" / 0.89 mm |
| Distal Tip: | Straight / 3mm J-Tip / 6mm J-Tip |
| Tip Configuration: | 5 cm / 10cm / 15cm Flexible Tip |
| Shaft Length: | 150cm / 180cm / 260 cm |

* a choice of various sizes will assist your patient's needs.

| LOKUM "Lunderquist" Order Codes (5 per Box) | | | |
|--|----------------|----------------|----------------|
| 0.035" / 0.89mm Tip / Shaft Length | Straight Tip | 3mm J-Tip | 6mm J-Tip |
| 5 cm / 150 cm | LOLU-05-150-ST | LOLU-05-150-3T | LOLU-05-150-6T |
| 5 cm / 180 cm | LOLU-05-180-ST | LOLU-05-180-3T | LOLU-05-180-6T |
| 5 cm / 260 cm | LOLU-05-260-ST | LOLU-05-260-3T | LOLU-05-260-6T |
| 10 cm / 150 cm | LOLU-10-150-ST | LOLU-10-150-3T | LOLU-10-150-6T |
| 10 cm / 180 cm | LOLU-10-180-ST | LOLU-10-180-3T | LOLU-10-180-6T |
| 10 cm / 260 cm | LOLU-10-260-ST | LOLU-10-260-3T | LOLU-10-260-6T |
| 15 cm / 150 cm | LOLU-15-150-ST | LOLU-15-150-3T | LOLU-15-150-6T |
| 15 cm / 180 cm | LOLU-15-180-ST | LOLU-15-180-3T | LOLU-15-180-6T |
| 15 cm / 260 cm | LOLU-15-260-ST | LOLU-15-260-3T | LOLU-15-260-6T |



AndraTec GmbH
Simmernerstr. 70
D-56075 Koblenz
Tel: +49 261 29675 172
Fax: +49 261 29675 173
www.andratec.com