



- The LOKUM Lunderquist **MINI** Guide Wire is especially designed for extra strong support and high pushability in difficult procedures.
- The high performance PTFE coating combined with extra stiff stainless steel LOKUM Lunderquist **MINI** is an excellent choice to treat difficult vascular and non-vascular lesions.
- The Lokum Lunderquist Type guide wire **MINI** is intended to use for complex interventional procedures where increased stability, and low surface friction of the guidewire is needed.

Technical Information:

Guide Wire Material:	Extra Tempered Stainless Steel
Coating:	High Performance PTFE
Diameter:	0.025" / 0.635 mm
Distal Tip:	Straight / 3mm J-Tip / 6mm J-Tip
Tip Configuration:	5 cm / 10cm / 15cm Flexible Tip
Shaft Length:	150cm / 260 cm

* a choice of various sizes will assist your patients needs.

LOKUM "Lunderquist-MINI" Order Codes (5 per Box)			
0.025" / 0.635mm Tip / Shaft Length	Straight Tip	3mm J-Tip	6mm J-Tip
5 cm / 150 cm	LUMI-05-150-ST	LUMI-05-150-3T	LUMI-05-150-6T
10 cm / 150 cm	LUMI-10-150-ST	LUMI-10-150-3T	LUMI-10-150-6T
15 cm / 150 cm	LUMI-15-150-ST	LUMI-15-150-3T	LUMI-15-150-6T
5 cm / 260 cm	LUMI-05-260-ST	LUMI-05-260-3T	LUMI-05-260-6T
10 cm / 260 cm	LUMI-10-260-ST	LUMI-10-260-3T	LUMI-10-260-6T
15 cm / 260 cm	LUMI-15-260-ST	LUMI-15-260-3T	LUMI-15-260-6T



AndraTec GmbH
 Simmernerstr. 70
 D-56075 Koblenz
 Tel: +49 261 29675 172
 Fax: +49 261 29675 173
 www.andratec.com