


# AndraTec *Advanced Interventional Technologies*



*PillaR*  
Bifurcation-Reconstruction-Device

- Unique Bifurcation Stent System to treat bifurcation narrowing or occlusion.
- Easy contralateral wire cannulation for contralateral stent insertion due to a unique entrance gate
- Alternative to kissing stents or multiple bifurcation stents
- Reconstruction Device intended to naturally re-establish the bloodflow in the vessels.
- Low profile delivery system with optimal wall coverage
- Advanced Remodeling

## Technical Information:

Material:	MP35N Cobalt Chromium
Diameter:	18 mm Recommended Max. Dilatation Range
Length:	61 mm
Profile:	PTA Balloon plus 1F / min 8F
Recommended Sheath:	8F
Recommended Guide Wire	0.035" (0.089 mm) Lokum Lunderquist AndraTec
Recommended Sheath Introducer:	Terumo destination

\* a choice of various sizes (diameters and lengths) will assist your patients needs.  
Future sizes will be available soon

Foreshortening Table PillaR					
Stent length (mm)	∅ 8	∅ 10	∅ 12	∅ 15	∅ 18
-	-	-	-	-	-
61	60	59	59	59	57
-	-	-	-	-	-

PillaR Bifurcation Reconstruction Device Order Codes			
Length Diameter		Length	Order Code
61 mm	Pillar Main Body (PMB)	61 mm	PR-25-61
28 mm	Pillar Aortic Extender (PE)	28 mm	PE-25-28
33 mm	Pillar Iliac Limb (PIL)	33 mm	PE-25-33

PillaR Covered B-R-D Order Codes			
Length Diameter		Length	Order Code
61 mm	Pillar 30mm PTFE covered	61 mm	PRCV-25-61

CE pending:

*Andratec*   
Advanced Interventional Technologies

AndraTec GmbH  
Simmernerstr. 70  
D-56075 Koblenz  
Tel: +49 261 29675 172  
Fax: +49 261 29675 173  
www.andratec.com

Rev. I / 2015-12