



Pillow
AAA-Occluding Device

- Long term results proof ; Unique Endovascular Aneurysm Repair and Thoracic Endovascular Aneurysm Repair occlusion Device
- Flexible and conformable within the Aneurysm sac
- Low profile delivery system
- Optimal Aneurysm sac filling to promote thrombosis and prevent ENDOLEAKS in variable sac locations
- The unique flexible and atraumatic design provides a barrier to blood flow thus enhancing intra-sac thrombosis
- Ease of use to ensure safe and atraumatic delivery and placement at the desired location
- Very low outward forces to prevent further pressure to the aneurysm sac

Technical Information:

| | |
|--------------------------------|-------------------|
| Material: | Braided Nitinol |
| Diameter: | 2.5 cm |
| Length: | 8 cm |
| Profile: | 8F |
| Recommended Sheath: | 9F / 90 cm |
| Recommended Guide Wire | 0.035" (0.089 mm) |
| Recommended Sheath Introducer: | Ranging 9F |

Pillow AAA Occlusion Device Order Codes

| Length Diameter | 40mm | 80mm | 120mm |
|--------------------|----------|----------|-------|
| 2.5 cm | PI-25-40 | PI-25-80 | - |

SDS Pillow Delivery System Order Code

| | |
|----------------------|----------|
| Length 160cm / 8F | SDS-1600 |
|----------------------|----------|

Rev. 1 / 2017



AndraTec GmbH
Simmernerstr. 70
D-56075 Koblenz
Tel: +49 261 29675 172
Fax: +49 261 29675 173
www.andratec.com