



lokum®
AT | LOKUM AMPLATZ TYPE
GUIDE WIRE

combining positioning and extra support in one design



- Design characteristics for positioning and extra support in one.
- High performance PTFE coated shaft providing low surface friction and enhanced trackability.
- Intended use during interventional peripheral and interventional biliary drainage procedures.
- Providing enhanced support during catheter and/or delivery system exchange.



Technical Information

Guide Wire Material	Extra Tempered Stainless Steel
Coating	High performance PTFE
Diameter	0.035" (0.89 mm)
Distal Tip	Straight / 3 mm J-Tip / 6 mm J-Tip
Tip Configuration	5 cm Flexible Tip
Shaft Length	180 cm / 260 cm / 360 cm

Order Codes

A choice of various diameters and lengths.

	3 mm 	6 mm 			
LOAM-180-ST	LOAM-180-3T	LOAM-180-6T	0.035" (0.89 mm)	180 cm	5 cm
LOAM-260-ST	LOAM-260-3T	LOAM-260-6T	0.035" (0.89 mm)	260 cm	5 cm
LOAM-360-ST	LOAM-360-3T	LOAM-360-6T	0.035" (0.89 mm)	360 cm	5 cm

Office & Warehouse

